



Wound Care

Compression Therapy for
Self-Care of Venous Ulcers



medi. I feel better.

mediven® dual layer stocking system

elastic, two-layer stocking system for venous ulcers.

- Dual layer stocking system provides accurate compression.
- Easy donning & doffing with closed toe liner and open toe stocking.
- Unisex knee length product.
- Firm & Extra Firm compression (30-40 & 40-50 mmHg).
- ClimaFresh odor protection inhibits bacterial growth within the liner.
- Under liner combines softness with ease of application.
- Outer stocking provides maximum wearing comfort to enhance compliance.
- Designed to work with a primary dressing.
- System includes: 2 closed-toe liners and 2 open-toe stockings.



Relevant Codes*:

HCPCS A6531 gradient compression stocking, below knee, 30-40 mmHg, each. Report with correct modifiers, RT (right) and/or LT (left).

HCPCS A6532 gradient compression stocking, below knee, 40-50 mmhg, each. Report with correct modifiers, RT (right) and/or LT (left).

Leg Circumference		SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
SIZE CHART	Calf	11" – 15" 28–38cm	12¼" – 16½" 31–42cm	13½" – 18" 34–46cm	14½" – 19½" 37–50cm	18" – 22" 46–56cm
	Ankle	7" – 8½" 18–21.5cm	8½" – 9¾" 21.5–25cm	9¾" – 11" 25–28cm	11" – 12½" 28–32cm	12½" – 14½" 32–37cm
Leg Length		Standard: 15"–19" (38–47.5cm)				

* Code usage should be verified by responsible parties in accordance with all published guidelines.

circaid® juxtalite® lower leg

inelastic, ready-to-wear compression garment for venous diseases.


- Inelastic, instantly adjustable compression.
- Multiple compression levels (20-30, 30-40, 40-50 mmHg).
- juxtalock™ band system allows for easy and instant adjustability of product without removal.
- Built-In Pressure system™ (BPS™) guide card provides targeted compression ranges for improved results.
- 6-month warranty.
- System includes: 1 legging, 1 pair of compression anklets, 1 undersleeve (footless), and 1 Built-In Pressure system (BPS) guide card.
- Does not include circaid power added compression band™ (pac band™) or juxtalite ankle foot wrap. Must be purchased separately.



Relevant Codes*:


HCPCS A6545 for an open venous stasis ulcer, billed with the AW modifier. Limited to one compression wrap per six months, per leg with RT (right) and/or LT (left) modifiers. Gradient compression, non-elastic, below knee 30-50 mmHg.

juxtalite lower leg



circumference in centimeters								
size	S	M	M - X	L	L - X	XL	XL - X	XXL
cC	26 - 36	33 - 44	44 - 54	41 - 51	51 - 61	37 - 48	48 - 61	53 - 64
cB	19 - 22	22 - 27	22 - 27	25 - 32	25 - 32	30 - 37	30 - 37	35 - 42
actual garment length in centimeters								
size	short				long			
lB - D	28				33			
sizing measurement			short			long		
l floor - D1			<44			>44		

All garments contain the circaid Built-In Pressure system (BPS) guide card.



circaid power added compression band (pac band)

circumferences in cm.			
	S	L	XL
cB	19 - 27	25 - 50	<56
cA1	19 - 36	22-41	<44

medi USA
6481 Franz Warner Parkway
Whitsett, NC 27377
mediusa.com
phone 1.800.633.6334
fax 1.888.570.4554