



Coverflex[®] grip tubular bandages

Coverflex[®] grip is an elastic tubular bandage, made from a cotton and synthetic blend, that does not contain natural rubber latex. It provides 360 degrees of uniform, consistent support and low level compression. The bandage is commonly used for secondary fixation, support for joints and ligaments, and to protect sensitive skin from adhesives or trauma.

Coverflex[®] grip sizing guide

Bandage Size	Bandage Dimensions (diameter x length)	Limb Measurements*			
		Min inches	Max inches	Min cm	Max cm
B	2.5" x 3 or 11 yds / 6.3cm x 2.7 or 10m	5.5	8	14	20
C	2.7" x 3 or 11 yds / 6.75cm x 2.7 or 10m	6	9	15	23
D	3" x 3 or 11 yds / 7.5cm x 2.7 or 10m	7	10	18	25
E	3.4" x 3 or 11 yds / 8.75cm x 2.7 or 10m	8	12	20	30
F	4" x 3 or 11 yds / 10cm x 2.7 or 10m	9	18	23	46
G	4.7" x 3 or 11 yds / 12cm x 2.7 or 10m	12	28	30	71
J	7" x 3 or 11 yds / 18cm x 2.7 or 10m	18	38	46	97
K	8.6" x 3 or 11 yds / 22cm x 2.7 or 10m	28	50	71	127



*measure largest circumference

Order Number	Bandage Size	Bandage Dimensions (diameter x length)	Selling UOM	Bandages per Box	HCPSC Code
931066	B	2.5" x 11 yds. / 6.25cm x 10m	Box	1	A6457
907000		2.5" x 3 yds. / 6.25cm x 2.74m	Bag		
931067	C	2.7" x 11 yds. / 6.75cm x 10m	Box	1	A6457
907001		2.7" x 3 yds. / 6.75cm x 2.74m	Bag		
931068	D	3" x 11 yds. / 7.5cm x 10m	Box	1	A6457
907002		3" x 3yds. / 7.5cm x 2.74m	Bag		
931069	E	3.4" x 11 yds. / 8.75cm x 10 m	Box	1	A6457
907003		3.4" x 3yds. / 8.75cm x 2.74 m	Bag		
931070	F	4" x 11 yds. / 10cm x 10m	Box	1	A6457
907004		4" x 3yds. / 10cm x 2.74m	Bag		
931071	G	4.7" x 11 yds. / 12cm x 10m	Box	1	A6457
907005		4.7" x 3yds. / 12cm x 2.74m	Bag		
931072	J	7" x 11 yds. / 18cm x 10m	Box	1	A6457
907006		7" x 3yds. / 18cm x 2.74m	Bag		
931073	K	8.6" x 11 yds. / 22cm x 10m	Box	1	A6457
907007		8.6" x 3yds. / 22cm x 2.74m	Bag		



All bandage lengths are unstretched

HARTMANN USA, Inc.
Rock Hill, SC 29730
1-800-243-2294

www.hartmannusa.com

For more information, contact your local sales representative
or call our Customer Service Center at 1-800-243-2294.

To learn more about HARTMANN products, please visit www.hartmannusa.com